

SurgiSeal Stylus®

TOPICAL SKIN ADHESIVE

2-Octyl Cyanoacrylate

PROFESSIONAL USE REFERENCE



1. Point the SURGISEAL Stylus® applicator towards the ceiling and away from the patient. Press the button of the applicator upward.



2. Invert the applicator and allow the adhesive to flow through the tip. No squeezing is necessary in normal application.

If you require greater flow, you may squeeze the ridged portion of the applicator.



3. Approximate the wound edges together and apply in one thick layer or two thin layers.

If two thin layers applied, allow 30 seconds dry time in between layers.



Medical Adhesive
Technologies

H.B. Fuller Medical Adhesive Technologies, LLC, 506 Pine Mountain Rd, Hudson, NC 28638, USA
Tel: +1.484.334.2929 · +1.866.579.4508 · www.hbfuller.com/medical
SurgiSeal Stylus® is a registered trademark of H.B. Fuller Medical Adhesive Technologies, LLC.
© 2025 H.B. Fuller Medical Adhesive Technologies, LLC. All rights reserved.

SSS-PUR1-2503